

COMPLETELY REUSABLE

A.T.T.S. SILO


Automatic Tissue Transfer System

1600 mL Collection

- Easy to assemble & disassemble
- Easy to clean & autoclave
- Easy to collect & re-inject
- Simple collection of harvested tissue
- Can be used “in-line” with aspirator
- Free Standing, or can be mounted to pole, or slide to your holder. (Optional)
- Disposable tubing with larger diameter also available
- Base & Lid made from Stainless steel
- Canister made from high temperature medical grade plastic and virtually indestructible



The base is tapered allowing tissue to draw smoothly



KMI IMI GROUP

4 Autry, Irvine,
CA, 92618, USA www.kigoc.com
Tel: +1(949) 458-1897

© 2018 all rights reserved. KMI IMI GROUP


26-0002

A.T.T.S. SILO

Automatic Tissue Transfer System

PARTS & ACCESSORIES

Item No.	Description
IM-ATTS-SILO-L	A.T.T.S. SILO Lid - Stainless Steel
IM-ATTS-SILO-C	A.T.T.S. SILO Canister - 1600cc
IM-ATTS-SILO-R	A.T.T.S. SILO Holding Rod - Stainless Steel (2 PCS)
IM-ATTS-SILO-B	A.T.T.S. SILO Base - Stainless Steel
IM-ATTS-SILO-BL	A.T.T.S. SILO Barbed Luer
IM-ATTS-SILO-LTL	A.T.T.S. SILO Luer Transfer
IM-ATTS-SILO-FLC	A.T.T.S. SILO Female Cap
IM-ATTS-SILO-12ST187	A.T.T.S. SILO 12" silicone Tubing, .187 Dia.
IM-ATTS-SILO-MLC	A.T.T.S. SILO Male Luer Cap
IM-ATTS-SILO-MB	A.T.T.S. SILO Mounting Bracket
IM-ATTS-SILO-SMB	A.T.T.S. SILO Slide Mounting Bracket
IM-ATTS-SILO-OB	A.T.T.S. SILO O-Ring for Base
IM-ATTS-SILO-OL	A.T.T.S. SILO O-Ring for Lid
IM-ATTS-SILO-BRU	A.T.T.S. SILO Cleaning Brush for Cleaning Base & Lid
IM-ATTS-SILO-YC	A.T.T.S. SILO Y-Connector - Stainless Steel
IM-ATTS-SILO-TC	A.T.T.S. SILO Plastic Tubing Clamp



KMI IMI GROUP

4 Autry, Irvine,
CA, 92618, USA www.kigoc.com
Tel: +1(949) 458-1897